ROUSELEY

ENG	Silicone tourniquet - Instruction For Use	2
FR	Garrot en silicone - Mode d'emploi	5
DE	Silikon-Tourniquet - Gebrauchsanweisung	8
IT	Laccio Emostatico in Silicone - Istruzioni per l'Uso	11
ES	Torniquete de Silicona - Instrucciones de Uso	14
PT	Torniquete de Silicone- Instruções de Utilização	17
PL	Opaska uciskowa Silikonowa (tourniquet) - Instrukcja użytkowania	20
NL	Siliconen Tourniquet - Gebruiksaanwijzing	23
SV	Silikon Tourniquet (Silikontorniké) - Bruksanvisning	26
NO	Silikon Tourniquet (Silikonturniké) - Bruksanvisning	29
FI	Silikoninen Torniketti (Kiristysside) - Käyttöohje	32
LT	Silikoninis Turniketas (Žgutas) - Naudojimo instrukcija	35
LV	Silikona Tournikets (Žņaugis) - Lietošanas instrukcija (IFU)	38
ET	Silikooni Turniquet (Silikoonist žgutt) - Kasutusjuhend	41
TR	Silikon Turnike - Kullanım Talimatı	44

Silicone tourniquet

Instruction For Use (IFU)

Product Name: Silicone Tourniquet
Model: TNQ-T0K00
Sizes: Adult and Pediatric
Classification: Class I, Non-Sterile, Reusable
Manufacturer: Rolseley Ltd, UK
EU Rep: Tikrola MB, Lithuania
Version & Date: v1.0 / 2025.06.25
Electronic Version: tourniquet.store/pages/IFU

ROLSELEY

Please carefully consider this instruction and read all the provisions prior to use. Follow the maintenance recommendations.

1. INTENDED USE & DESCRIPTION

Purpose: Temporary restriction of venous flow during blood draws or IV access.

Indications: Venipuncture, blood collection, IV insertion.

Contraindications: Circulatory disorders, silicone allergy, damaged skin, prolonged application.

Target Users: Trained healthcare professionals only.

Patient Groups: Adults and pediatric patients (select correct size).

Use Environment: Hospitals, clinics, labs, ambulances.

Functional Principle: The tourniquet works by applying controlled pressure to a limb using an elastic silicone band and buckle mechanism to temporarily restrict venous blood flow.

Expected Clinical Benefits: Enhances vein visibility for safer venipuncture, reduces needlestick attempts, supports routine clinical use, minimizes skin reactions (latex-free), and maintains reliable performance with proper care.

Size: Available in Adult and Pediatric sizes. Select appropriate size based on patient limb circumference.

Dimensions: Adult size: 58 cm x 2,5 cm; Pediatric size: 42 cm x 1,5 cm.

Accessories & Spare Parts: None required.

Device Compatibility: Not to be used with automated or electronic compression devices.

2. APPLICATION & USE

Preparation for Use:

This device requires no installation, assembly, or calibration.


Before use, check for:

- Cracks, cuts, loss of elasticity.
- Buckle malfunction.
- Residue that cannot be removed.

△ Discard if damaged.

Application

Place 5–7 cm above puncture site on clean, intact skin (fig. 1). Tighten to restrict venous, not arterial, flow (fig. 2). Avoid pinching or over-tightening. Remove immediately after use by pressing the release button (fig. 3).

Warning: Never apply on or near the neck — risk of suffocation. Any use beyond the stated intended purpose is improper and at the operator's own risk. The manufacturer is not liable for such misuse.

3. CLEANING & REUSE

Before First Use:

- Disinfect with alcohol based wipes.

After Each Use:

- Clean manually with mild soap & water.
- Disinfect all surfaces with alcohol based wipes.
- Air dry completely.

Do Not:

- Autoclave, boil, or microwave.
- Exceed 100 °C during cleaning or drying.
- Use abrasives, bleach, or high-pH disinfectants.
- Use automated or machine-based cleaning systems.

Important:

Do not use disinfectant wipes that contain strong oxidizing agents, quaternary ammonium compounds, or high pH levels. Some commercially available wipes (e.g., **Clinell Universal Wipes**) have been reported by users to degrade the silicone surface over time and are not recommended.

4. STORAGE & SHELF LIFE

Shelf life: 3 years.

Storage: Dry, cool place (10–30°C), away from sunlight and chemical exposure.

Recommended lifespan: up to 1,500 uses or 3-6 months from first use.

5. RISK & SAFETY INFORMATION

Improper use may cause bruising, numbness, or nerve damage. Temporary skin marks or discomfort may occur; discontinue use if pain, redness, or swelling persists. Risks are deemed acceptable (AC/AFARP) per Risk Management Report CE/MDR-TNQ2025-07 / UK/TNQ2025-07.

⚠ Warnings & Precautions:

- Do not apply for over 1 minute. If unsuccessful, wait 2 minutes before retrying.
- Monitor skin color, temperature, and compression effect.
- Use only on limbs—not near or on the neck (risk of suffocation).
- Cut immediately if stuck and removal is unsafe.
- Disinfect after each use to prevent bacterial contamination.
- Avoid use if damaged, discolored, or altered; dispose of properly.
- Protect from heat, sunlight, and chemicals.
- Not for use by untrained persons, children, or impaired individuals.
- Rare silicone allergies may occur—discontinue if irritation appears.

6. DISPOSAL

General waste if disinfected.

Medical waste if contaminated.

Dispose of per local laws and EU Waste Directive 2008/98/EC.

7. PACKAGE & LABEL SYMBOLS

Contents: 1x Silicone Tourniquet.

Packaging: Sealed pouch inside paper box.

Sterility: Non-sterile; disinfect before first use.

Storage: See Section 4.

UDI: The Unique Device Identifier (UDI) specific to this product variant (e.g., size and color) is printed on the product label and packaging.

Labeling Symbols (CE & UKCA marked; compliant with EN ISO 15223-1:2021):

	CE conformity mark		UKCA conformity mark
	Reference number		Medical device
	Manufacturer		Lot code
	Attention		Date of manufacture
	Keep away from moisture		Please follow Instruction for use
	EU representative		Keep out of sunlight
	Do not use if package is damaged		Use-by date
	Multiple patients multiple use		Non sterile product
	Temperature limit - 10°C / 30°C		Unique Device Identifier

8. CONTACT

Manufacturer information:


ROLSELEY LTD
38 Arundel street,
Brighton, BN2 5TH
United Kingdom
+44 7584 655333
info@rolseley.co.uk

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Incident Reporting:

Report any serious incident or field safety corrective action to the Manufacturer (Rolseley Ltd) and, as applicable, to the MHRA (UK) or the national competent authority of the EU Member State where the incident occurred, in accordance with UK MDR 2002 and EU MDR 2017/745 Article 87.

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Yerel yasalara ve **AB Atık Direktifi 2008/98/EC** hükümlerine uygun şekilde imha ediniz.
AB Temsilcisi: